



Lemon-E

Detergent/Disinfectant

Detergente/Desinfectante

- Concentrated • Lemon Fragrance
- Cleans, Deodorizes & Disinfects Hard Surfaces
- EPA Registered for Sale In New York & New Jersey



Technical Data Sheet

PRODUCT DESCRIPTION:

An economical lemon scented quat disinfectant that has a broad spectrum of claims against gram negative and gram positive bacteria. It is especially effective at controlling odors that are bacterial in nature.

USES:

Use on inanimate surfaces such as glass, plastic, metal, porcelain, tile, floors, walls, appliances and furnishings. Use in restrooms, hospitals, schools, public rooms, homes, pet quarters, industrial and commercial offices and areas.

ADVANTAGES:

Pleasant lemon fragrance, cost effective, disinfects, cleans, deodorizes, non-corrosive, non-staining to plastic, vinyl, tile, enamel, synthetics and most common metals. Excellent for use in nursing homes where a low odor disinfectant is desirable.

SPECIFICATIONS:

Active Ingredients: 1.562%
 pH: 12
 Color: Yellow
 Solubility in Water: Complete
 Odor: Lemon
 Flash Point: None
 Viscosity: Water Thin
 Weight per Gallon: 8.33 Lbs.

PRODUCT/CONTAINER

DISPOSAL:

Dispose of used/unwanted product according to state and local regulations. Offer container for recycling.

SIZES AVAILABLE:

1 Gallon Bottle, 5 Gallon Pail, 30 & 55 Gallon Drums.



Recyclable packaging

DIRECTIONS FOR PRODUCT USE:

DISINFECTING: Before disinfecting, pre-clean surfaces as directed on product label. Then, prepare 4 oz. of Lemon-E per gallon of water and thoroughly wet the surfaces to be disinfected. The disinfecting solution must contact the surface for 10 minutes for proper disinfection. If desired, the surfaces may then be rinsed with clear water or wiped with a damp cloth.


CLEANING AND DEODORIZING: Use 4 oz. of Lemon-E per gallon of water. Cleaning may be accomplished with a clean mop, sponge, squeegee, wet vacuum system or a mechanical spraying device. After cleaning, remove suspended soils with a clear water rinse, vacuum pickup or by wiping with a clean mop, sponge, squeegee or cloth.

WHEN USED AT THE DILUTIONS SPECIFIED ON THE PRODUCT LABEL, LEMON-E IS EFFECTIVE AGAINST:

Escherichia coli, Staphylococcus aureus Methicillin resistant, Klebsiella pneumoniae, Legionella pneumophila, Staphylococcus aureus Vancomycin resistant — VISA, Pseudomonas aeruginosa, Salmonella enterica, Staphylococcus aureus, Salmonella schottmuelleri, Streptococcus salivarius, Adenovirus Type 2, Herpes simplex Type 2, Herpes simplex Type 1, HIV-1, Human Coronavirus, Avian Influenza A virus, Influenza A H1N1 virus (formerly known as Swine flu), Influenza A2/Japan, Hepatitis B Virus (HBV), Hepatitis C Virus (HCV), Trichophyton mentagrophytes, Aspergillus niger, Staphylococcus aureus, Klebsiella pneumoniae.



It is recommended that protective gloves and safety glasses be worn when using this and all cleaning products. Refer to MSDS for further safety information.



MANUFACTURED BY:



39-30 Review Avenue
 Long Island City, New York 11101
 1-888 WAX PLUS • (718) 392-8080 Fax: (718) 482-9366
 Visit us on the web! www.cleaning-solutions.com

